

Veritas microdissection system: optimized protocol for laser microdissection of living *in vitro* cells



VERITAS APPLICATION NOTE

By Steven Blakely and Shirley Chu, Molecular Devices Corporation, 1311 Orleans Dr. Sunnyvale, CA 94086.

INTRODUCTION

Laser Capture Microdissection (LCM) is a proven technique for the isolation of pure cell populations for downstream molecular analysis. The combined use of UV laser cutting with LCM using the infrared (IR) laser allows for the rapid and precise isolation of larger numbers of cells, while maintaining cellular and nucleic acid integrity necessary for downstream analysis. In this note, it is shown that these established techniques can also be used for the isolation of living cells, avoiding other more laborious methods of cell selection, and providing opportunity for a wide range of research applications.

This application note describes a simple and effective protocol for the isolation of living adherent cells and the successful subsequent re-cultivation of homogeneous subpopulations. To demonstrate the success of this technique, two cell lines were used: TM3 (normal mouse testis) and SKBR3 (human breast cancer).

MATERIALS

- Arcturus Veritas™ Microdissection Instrument (Molecular Devices Cat. #VERITAS)
- PEN membrane frame slide (Molecular Devices Cat. #LCM0521)
- CapSure® Macro LCM caps (Molecular Devices Cat. #LCM0211 or #LCM0212)
- Trypsin-EDTA (0.25% Trypsin, 1 mM EDTA, Invitrogen Cat. #25200-056)
- Hanks' Balanced Salt Solution (Invitrogen Cat. #14170-161)
- Petri dish, sterile, 100 mm diameter (Falcon Cat. #351005)
- Cover glass, 24 x 60 mm, No. 1 thickness (VWR Cat. #48393-106)
- RNase Away® (Life Technologies Cat. #10328-011)
- 100% ethanol (VWR Cat. #34172-00)

- CO₂ incubator (MLS)
- Pipettor and sterile pipette tips (MLS)
- 2-well chamber slides (VWR Cat. #62407-325)
- ND-1000 UV/Vis Spectrophotometer (NanoDrop Technologies)
- Agilent 2100 Bioanalyzer (Agilent Technologies Cat. #G2940CA)
- SKBR3 cell culture medium
 - McCoy's modified 5A media (Invitrogen Cat. #16600-082)
 - 10% Fetal Bovine Serum (FBS, Invitrogen Cat. #16600-044)
 - 1 mM Sodium Pyruvate (Invitrogen Cat. #11360-070)
 - 10 U Penicillin/10 µg Streptomycin (BioWhittaker Cat. #17-60ZE)
- TM3 cell culture medium
 - Dulbecco's Modified Eagle's Medium (Invitrogen Cat. #30-2002)/Ham's F12 Medium (Invitrogen Cat. #117650-054), 1:1
 - 5% Horse Serum (Invitrogen Cat. #16050-122)
 - 2.5% Fetal Bovine Serum (Invitrogen Cat. #16600-044)
 - 10 U Penicillin/10 µg Streptomycin (BioWhittaker Cat. #17-60ZE)

METHODS

Specimen preparation

Step 1. Thoroughly rinse the metal frame PEN membrane slide with 100% ethanol and air-dry prior to use. Keep slide in a sterile environment.

Note: Ensure the slide is completely dry prior to use.

Step 2. Trypsinize adherent cells from growth vessel (plate, flask) using standard protocol.

Step 3. Deactivate trypsin with media using standard protocol.

Step 4. Resuspend 1–2 mL of trypsinized cells in 10 mL of fresh media.

Step 5. Place a metal frame membrane slide with chamber facing up into a sterile Petri dish. Transfer 1 mL of the cell suspension into the chamber of the frame membrane slide. If necessary, gently rock the slide in the Petri dish to completely cover the chamber area with media.

Step 6. Place lid onto the Petri dish and incubate using appropriate culturing conditions for the cells until desired cell confluency is achieved.

Note: Replace with fresh media as needed.

Laser microdissection slide preparation

Step 1. Thoroughly clean the instrument and work area, including pipettors, pipette tip box, etc. with 100% ethanol and RNase-Away or RNase-Zap®.

Step 2. Thoroughly rinse a cover glass with 100% ethanol and air-dry prior to use. Keep cover glass in a sterile environment.

Note: Ensure the cover glass is completely dry prior to use.

Step 3. When cells have reached the desired confluency, remove the media from the chamber using a clean pipette tip.

Step 4. Add 950–1000 µL of fresh media to the chamber.

Step 5. Carefully place a cover glass over the chamber side of the frame slide to create a mini-environment for the cell culture, enabling extended survival and reducing the possibility of the cells drying out. (See Figure 1.)

Note: Take care to reduce the amount of air bubbles formed when applying the cover slip.

Step 6. Using a Kimwipe®, carefully blot any excess media that has seeped outside the cover glass.

Step 7. Transport slide inside the Petri dish to the Veritas system.

Step 8. Remove the slide from the Petri dish and use a Kimwipe soaked in 100% ethanol to clean the flat side of the frame slide. Dry the slide completely.

Note: Make sure not to rupture the membrane.

Step 9. Insert the frame slide with the chamber and cover glass facing down (flat side up) onto the Veritas instrument and proceed to the laser microdissection session. (See Figure 2.)

Laser microdissection protocol

The following settings were used for protocol validation and should be used as a guideline for the microdissection of live cells. Optimization of settings may be required depending on the individual cell preparation.

Important:

→ Use CapSure Macro LCM caps.

→ Perform cut and capture at 10x or 20x.

→ Always capture first and then cut.

→ Turn the visualizer off.

Cutting (UV) laser settings:

→ UV Laser Power: 20–25

→ UV Spacing: 5000 µm

→ Tab size: 1

→ Automatic LCM spots: None


Capture (IR) laser settings:

→ IR laser power: 80 mW

→ Pulse: 4000 ms


→ LCM spot overlap: 40%

laser microdissection cell culture chamber (figure 1)



Cells cultured in “well” of a metal frame membrane slide. The live cell culture is covered by a large cover glass prior to microdissection.

inverted laser microdissection cell culture chamber with CapSure LCM cap (figure 2)



The CapSure cap sits atop the inverted chamber in contact only with the PEN membrane during microdissection.

Step 1. Locate cells of interest to capture.

Step 2. Use the Cut Line feature to draw around cells.

Step 3. Use the Single Point capture feature to apply LCM spots that will fuse LCM membrane to PEN membrane.

Note: Apply adequate amount of LCM spots for the given region. Try to avoid placing spots directly onto areas containing cells.

Step 4. Place the CapSure Macro LCM cap onto the area of the slide containing cells of interest.

Step 5. Locate LCM laser and fire a test LCM shot. If necessary, adjust laser settings.

Note: Confirm that LCM film has made contact to PEN film. LCM spot will be dark.

Step 6. Locate UV cutting laser.

Step 7. Activate LCM laser first and then the UV cutting laser.

Step 8. Move Macro LCM cap to QC station.

Step 9. Confirm presence of cells on LCM cap. (See Figure 3.)

Step 10. Move cap to offload station.

Reculturing of captured live cells

Step 1. Remove Macro LCM cap from the offload station and invert. Place cap with isolated cells facing up into a clean Petri dish.

Step 2. Pipette 50 μ L of Hank's solution onto the Macro LCM cap film surface. Pipette up and down 2–3 times, and dispose of the solution.

Step 3. Pipette 50 μ L of Trypsin-EDTA directly onto the captured cells on the cap and incubate for 5 minutes at room temperature. Cover with Petri dish lid during this incubation.

Step 4. After incubation, pipette the Trypsin-EDTA up and down several times to ensure a single cell suspension, then transfer the cell suspension into a well of a sterile chamber slide (or alternate desired growth vessel) containing 1–2 mL of appropriate cell media.

Step 5. Place the chamber slide in the incubator under appropriate conditions. Wait 2–3 days for detection of cell growth. Monitor cell growth using standard culture technique, changing media as needed.

Step 6. Recultured cells may be used as desired for further experiments.


RESULTS AND CONCLUSIONS

Using two different cell lines, TM3 and SKBR3, crop circles were captured using the above techniques, each onto a separate CapSure LCM cap. Captured areas were designated as “Large” (2–3 mm diameter) or “Small” (500–700 μ m diameter) and contained varied numbers of cells. Subsequent to laser microdissection, some caps containing live cells were recultured following the above steps for protocol validation.

Alternately, microdissected cells were lysed and RNA was extracted to evaluate if there was any negative effect on the nucleic acid as a result of the microdissection protocol. These caps, containing captured live cells, were immediately placed into tubes containing PicoPure® RNA extraction buffer. The caps were placed into the extraction buffer prior to any reculture steps, immediately following removal of the cap from the Veritas microdissection instrument. Total RNA was isolated following the standard PicoPure protocol, and eluted in 30 μ L of elution buffer.

Total RNA yield for each sample was assessed using the NanoDrop ND-1000. (See Table 1.)

microdissection of living *in vitro* cells using the Veritas microdissection instrument (figure 3)



A: Live cells in media before laser microdissection. **B:** Cell culture chamber after microdissection. **C:** Microdissected live cells on the CapSure Macro LCM cap.

The total RNA yields varied between replicates, as the cell counts within the captured areas were not monitored. However, the larger crop circles did result in higher yields than the smaller areas, as expected, and the RNA yield for each captured area, regardless of size, yielded sufficient RNA for downstream molecular analysis. The amount of RNA obtained from live cells subsequent to laser microdissection also indicates that the isolation protocol did not have a negative effect on the nucleic acids.


RNA profiles were obtained using the Agilent 2100 Bioanalyzer. (See Figure 4.) Clear 18S and 28S peaks are seen in each profile, further confirmation that the microdissection process does not affect RNA integrity.

Table 1. Total RNA yield from microdissected live cells

Sample ID	Cell Type	RNA Conc. ng/ μ L	Total RNA Yield (ng)	Average Yield (ng)
Large Crop 1	SKBR-3	38.73	1161.90	964.65
Large Crop 2	SKBR-3	25.58	767.40	
Small Crop 1	SKBR-3	4.35	130.50	131.64
Small Crop 2	SKBR-3	6.73	201.90	
Small Crop 3	SKBR-3	3.77	113.10	
Small Crop 4	SKBR-3	4.44	133.20	
Small Crop 5	SKBR-3	2.65	79.50	
Small Crop 1	TM-3	4.34	130.20	184.74
Small Crop 2	TM-3	1.13	33.90	
Small Crop 3	TM-3	7.78	233.40	
Small Crop 4	TM-3	7.14	214.20	
Small Crop 5	TM-3	10.40	312.00	

Note: 1.5 μ L of each total RNA sample was measured using the NanoDrop ND-1000 UV/Vis spectrophotometer.

Bioanalyzer profiles of total RNA isolated from microdissected live cells (figure 4)



An amount of 1 μ L of each sample was loaded onto the Bioanalyzer and run per the manufacturer's protocol. Note: TM3 Small Crop #5 is not present due to sample loading error.

SALES OFFICES

United States

Molecular Devices Corp.
Tel. +1-800-635-5577
Fax +1-408-747-3601

Australia

Molecular Devices Pty. Ltd.
Tel. +61-3-9896-4700
Fax +61-3-8640-0742

Brazil

Molecular Devices Brazil
Tel. +55-11-3616-6607
Fax +55-11-3871-9994

China

Molecular Devices Beijing
Tel. +86-10-6410-8629
Fax +86-10-6410-8601

Molecular Devices Shanghai

Tel. +86-21-6887-8820
Fax +86-21-6887-8890

Germany

Molecular Devices GmbH
Tel. +49-89-9605-880
Fax +49-89-9620-2345

Japan

Nihon Molecular Devices Osaka
Tel. +81-6-6399-8211
Fax +81-6-6399-8212

Nihon Molecular Devices Tokyo

Tel. +81-3-5282-5261
Fax +81-3-5282-5262

Norway

Molecular Devices Skatron
Tel. +47-32-859-800
Fax +47-32-859-801

South Korea

Molecular Devices Korea, LLC
Tel. +82-2-3471-9531
Fax +82-2-3471-9532

United Kingdom

Molecular Devices Ltd.
Tel. +44-118-944-8000
Fax +44-118-944-8001

www.moleculardevices.com

CapSure and PicoPure are registered trademarks and Veritas is a trademark of Molecular Devices Corporation. All other trademarks are the property of their respective owners.
©2006 Molecular Devices Corporation. Printed in U.S.A.
7/06 EP #0120-1423A

